IEM®

HMS CS

DA

DE

EL

EN

ES

FR

HR

IT

NL

PL

PT

SV

TR

Brugsanvisning

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⁷ Fantin F. et al. Is augmentation index a good measure of vascular stiffness in the elderly? Age and Ageing 2007; 36: 43-48.

⁸ Nunan D. et al. Assessment of central haemodynamics from a brachial cuff in a community setting. BioMed Central 2012; 1471-2261.

⁹ The Task Force for the Management of Arterial Hypertension of the European Society of Hypertension (ESH) and the European Society of Cardiology (ESC). 2007 Guidelines for the management of arterial hypertension. European Heart Journal 2007; 28: 1462-1536.

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Instructions for Use
EN

Instructions for Use

Hypertension Management Software Client Server (HMS CS)

Version 6.2 and above

HMS CS is used to analyze measurements taken using the Mobil-O-Graph® or Tel-O-Graph®.

Caution: Federal law restricts this device to sale by or on the order of a physician

IEM GmbH
Gewerbepark Brand 42
52078 Aachen
Germany

E-Mail: info@iem.de
Website: www.iem.de

The contents of these instructions for use must not be duplicated or published without the written authorization of IEM GmbH.

The analysis software is protected by copyright law and is the property of the manufacturer. All rights reserved. The analysis software must not be uploaded, copied, decompiled, redeveloped, disassembled or brought into any human-readable format. All usage and ownership rights to the software remain with IEM GmbH.

© IEM GmbH 2023. All Rights reserved.

Version 5.1 – 13.03.2023 – EN

Table of contents

1	Introduction.....	7
1.1	Intended Use	7
1.2	Improper use.....	8
2	Safety and security.....	9
2.1	Definition of safety symbols.....	9
2.2	Important safety and security information.....	10
3	Description of HMS CS.....	11
4	Working with the HMS CS.....	12
5	Installing the Software	13
5.1	System requirements	14
5.2	Installing on Windows®	15
5.2.1	Installing HMS CS from website.....	15
5.3	Installing on macOS.....	16
5.3.1	Installing HMS CS from website.....	16
5.4	Updating HMS CS.....	16
6	Toolbar.....	17
7	Starting and exiting HMS CS	18
8	First steps with the sample patient.....	19
8.1	Patient information	19
8.2	Blood pressure–PWA	20
9	Editing patient information.....	21
9.1	Signing-in a new patient.....	21
9.2	Selecting an existing patient	21
9.3	Adding to and changing patient information	22
9.4	Defining individual blood pressure limits	22
9.5	Deleting a patient.....	22
9.6	Importing blood pressure readings manually	23
10	Program configurations (options)	24
10.1	General	24
10.1.1	General	24
10.1.2	Units.....	24
10.1.3	Calibration.....	25
10.1.4	Language	25
10.1.5	Database.....	25
10.1.6	Global blood pressure limits	26
10.1.7	About	27

10.1.8	Import ABDM.mdb	28
10.1.9	Import GDT	28
10.1.10	Import patients	28
10.1.11	Data backup	28
10.1.12	Restore data.....	29
10.2	Port settings	29
10.2.1	Bluetooth® interface.....	30
10.2.2	Serial/USB interface.....	30
10.2.3	Deleting the monitor from the list	31
10.3	Report	31
10.4	GDT settings.....	32
10.5	Export.....	33
10.5.1	Export file name.....	33
10.5.2	Formatting	33
10.5.3	Exporting as CSV.....	33
10.6	Audit-Trail.....	33
10.7	Customization.....	33
10.7.1	Printing	33
10.7.2	Logos.....	34
10.7.3	Colors	34
11	Printing.....	35
11.1	Printing in-office blood pressure measurement.....	35
11.2	Printing 24-hour blood pressure measurement	35
12	Exporting measurement data	36
12.1	Exporting a patient's complete data	36
12.2	Exporting individual measurements	36
13	Data analysis.....	37
13.1	Measurement data analysis	38
13.2	Entering findings for measurement series	38
13.3	Blood pressure – PWA tab	39
13.3.1	Table of measurements.....	39
13.3.2	Blood pressure profile	40
13.3.3	Bar chart.....	42
13.3.4	Report	43
13.3.5	Hourly mean values	43
13.4	Comparing multiple measurement series	44
13.4.1	24h PWA	44

13.4.2	Patient report	45
13.4.3	Arterial stiffness	46
13.5	ECG tab.....	47
13.6	Weight tab.....	48
14	HMS CS & Mobil-O-Graph®	49
14.1	Bluetooth® connection	49
14.1.1	Pairing the blood pressure monitor with HMS CS	49
14.1.2	Prepare the blood pressure monitor for long-term measurements	51
14.1.3	Start the 24-hour measurement.....	52
14.1.4	Transferring and storing long-term measurement results.....	52
14.2	Connection via cable.....	54
14.2.1	Connecting the blood pressure monitor to the computer via cable.....	54
14.2.2	Configure the blood pressure monitor in HMS CS.....	55
14.2.3	Prepare the blood pressure monitor for long-term measurements	56
14.2.4	Start the 24-hour measurement.....	57
14.2.5	Connect the blood pressure monitor to the computer again	58
14.2.6	Transferring and storing long-term measurement results.....	59
15	HMS CS & Tel-O-Graph® BT.....	60
15.1	Connection via Bluetooth®	60
15.1.1	Pairing the blood pressure measuring device with HMS CS	60
15.1.2	Prepare the blood pressure measuring device for blood pressure measurement	62
15.1.3	Start blood pressure measurement.....	63
15.1.4	Transferring and storing blood pressure measurement results	64
16	Prepare device	66
16.1	Patient list	66
16.2	Setting device clock	66
16.3	Importing patient ID	66
16.4	Testing the device	67
16.5	Deleting old measurements/records.....	67
16.6	Special Mobil-O-Graph® settings	68
16.6.1	Configuring the log.....	68
16.6.2	Pre-set	68
16.6.3	In-office blood pressure monitoring	68
17	Pulse wave analysis	71
17.1	In-office pulse wave analysis	72
17.1.1	PWA license key for Mobil-O-Graph® / Tel-O-Graph® BT	72
17.2	PWA license dongle for Mobil-O-Graph® (Not available in the USA).....	73

17.3 In-office pulse wave analysis	73
17.3.1 Individual pulse wave analysis.....	73
17.3.2 Triple pulse wave analysis.....	74
17.4 24-hour PWA	76
17.4.1 Performing 24-hour PWA using Mobil-O-Graph [®]	76
17.4.2 Uploading and analysis of the 24-hour PWA.....	77
17.5 Displaying the pulse wave analysis.....	77
17.5.1 Alx and PWV history	79
18 Error messages	80

1 Introduction

Thank you for choosing Hypertension Management Software Client Server (HMS CS). These instructions for use will quickly familiarize you with the use of HMS CS and the many applications it can be used for.

HMS CS allows

- measurements to be transferred via a serial/USB interface or via Bluetooth®
- measurements to be evaluated for analysis purposes using various presentational methods and diagram types
- patient data to be managed.

The following IEM® products can be used in combination with HMS CS:

- Mobil-O-Graph® 24-hour blood pressure monitor
- Tel-O-Graph® blood pressure measuring device

► Note

The ECG functionality is no longer supported in this HMS CS version.

Therefore, the connection of the BEAM® with the HMS CS 6.2.0 and newer is no longer possible.

Furthermore, the ECG and weight data visualizations are no longer available.

Further notes, also on retrieving old ECG and weight data, can be found in Chapter 13.5 ECG tab and 13.6 Weight tab.

Read these instructions for use carefully before use and keep it in a suitable place so that the information is always available to you when you require it.

► Note

If you require a printed copy of the manual, please contact your specialist supplier or IEM® GmbH.

1.1 Intended Use

The HMS CS software and related accessories are used to prepare medical equipment for the purpose of measuring the blood pressure and vascular state of a patient. The software is intended to be used to measure, analyze, format, display, print and store cardiovascular parameters to assist in and monitor diagnoses and hypertension. The analysis of the blood pressure parameters allows the doctor to adjust the patient's medication accordingly and provide the patient with relevant lifestyle advice. HMS CS may only be used by or following the instructions of a medical doctor. Only a medical doctor may analyze the recorded and displayed HMS CS data.

The following medical devices may be used with the HMS CS software in keeping with the device's "intended use":

- Combined with the Mobil-O-Graph® blood pressure monitor, the "intended use" would be 24-hour blood pressure monitoring and pulse wave analysis (PWA).
- The Tel-O-Graph® blood pressure measuring device can be used with HMS CS and is used for blood pressure monitoring and pulse wave analysis (PWA).

For Tel-O-Graph[®] devices Pulse Wave Analysis (PWA) is not available in the USA. In general the following parameters are not available in the USA:

- Pulse pressure amplification
- Stroke volume (SV)
- Cardiac output (CO)
- Total vascular resistance (TVR)
- Cardiac index
- Augmentation pressure
- Pulse wave velocity (PWV)

WARNING

Please refer to the relevant instructions for use for further important information on the medical devices Mobil-O-Graph[®] and Tel-O-Graph[®] and their functions!

1.2 Improper use

HMS CS must not be used for any purpose other than the blood pressure analysis and data management procedures described here.

HMS CS in conjunction with Mobil-O-Graph[®] or Tel-O-Graph[®] is not designed to provide monitoring services with alarm signals in intensive care units or elsewhere.

There is no proven scientific knowledge about the use of HMS CS with new-born infants, pregnant women or for pre-eclampsia.

WARNING

Please refer to the relevant instructions for use for further important information on individual products (Mobil-O-Graph[®] and Tel-O-Graph[®]) and their functions!

2 Safety and security

Read the safety information carefully before using the products! It is important that you understand the information in these instructions for use. Please contact Technical Support if you have any questions.

2.1 Definition of safety symbols

The following symbols and signal words are used in this manual to indicate dangers and important information:

WARNING

Short description of the danger

This warning symbol in connection with the signal word **WARNING** indicates a possible or immediately threatening danger.

Failure to comply may lead to minor, moderate or severe injuries or death.

CAUTION

Short description of the danger

This warning symbol, in connection with the signal word **ATTENTION**, indicates possible material damage.

Non-adherence may lead to damage to the products or their accessories.

Note

The **NOTE** signal word indicates further information on HMS CS.

External Reference

This symbol indicates references to external documents that provide additional optional information.

Tip

This symbol refers to useful tips that may help you work.

2.2 Important safety and security information

WARNING

This system must not be used to provide blood pressure monitoring services with alarm signals in operating theatres or intensive care units.

CAUTION

To ensure data security, please observe the following:

- Do not set up a guest account on the computer.
- Use the HMS CS data backup function for regular backups. HMS CS does not create automatic backups.
- Regularly update your operating system, firewall and anti-virus software.
- Do not use operating systems for which support has been discontinued.
- Ensure that only authorized personnel have access to your computer.

CAUTION

Please take note of the battery charge level in the blood pressure measuring device when performing long-term blood pressure measurements. If the charge is not sufficient, a relevant message will appear in HMS CS!

Note

The use of HMS CS requires a basic knowledge of and experience with the operating systems Windows® and macOS.

Note

- The pulse wave analysis provides additional indicators of potential risks, but it is not a sufficiently reliable indicator of individual illnesses or treatment recommendations.
- Please note that the use of pulse wave analysis on children is not currently supported by clinical studies using reference methods.

3 Description of HMS CS

Blood pressure is measured, exported and stored in HMS CS, where you can analyze the measured values in line with your needs.

The patient information includes all key data such as:

1. Patient ID (mandatory input)
2. Name (mandatory input)
3. Address, telephone
4. Personal data (age, gender, etc.)
5. Medication, medical history, emergency contacts

HMS CS offers you various options to analyze the 24-hour ABPM. Results can be displayed on your computer screen, selected or printed out:

- Display all individual measurements
- Statistical analysis with mean blood pressure values for the entire day, daytime and night-time, as well as mean hourly values
- Extreme values (maximum, minimum)
- Graphical analyzes:
 - Envelope curve of mean hourly values
 - Bar chart of measurements
 - Curve of measurement values
 - Curve comparison for treatment optimization

4 Working with the HMS CS

Note

The use of HMS CS requires a basic knowledge of and experience with the operating systems Windows[®] and macOS.

HMS CS is used to manage and analyze blood pressure measurement data. These measurement values are assigned to a patient. Any number of measurement series can be stored for any given patient.

The following steps are performed in HMS CS:

Prepare measurement

- Start HMS CS
- Select existing patient or sign-in new patient
- Connect monitor to HMS CS
- Prepare device

Process measurement data

- Start HMS CS
- Select patient
- Connect HMS CS to monitor
- Read-out measured values from device
- Analyze measurement data
- Close HMS CS.

5 Installing the Software

The HMS CS software can communicate with the monitor over a variety of interfaces. The connection options are:

- Bluetooth®,
- Cable with serial interface (e.g. COM1, COM2),
- Cable with USB interface,

Note

Bluetooth® is not supported on operating system macOS.

Re-installation:

Download the HMS CS from the following URL: www.iem.de/hmscs.

Software update:

We recommend performing HMS CS software updates through our local trained partners.

Note

When using a USB cable, please install the USB driver before connecting the cable to the computer.

CAUTION

To ensure data security, please observe the following:

- Do not set up a guest account on the computer.
- Use the HMS CS data backup function for regular backups. HMS CS does not create automatic backups.
- Regularly update your operating system, firewall and anti-virus software.
- Do not use operating systems for which support has been discontinued.
- Ensure that only authorized personnel have access to your computer.

5.1 System requirements

Computer

- Processor: min. 1 GHz
- Memory: min. 2 GB RAM
- Hard disk space: min. 500 MB
- Resolution: min. 1024x768 pixels
- At least one free USB socket

Operating System

- Windows[®] 10 (32-bit & 64-bit), Windows[®] 11
- macOS 13 Ventura (without Bluetooth[®] support)

Note

Bluetooth[®] is not supported on operating system macOS.

Software

- Java Runtime Environment

Bluetooth[®]

- Bluetooth[®] USB adapter
- Bluetooth[®] 2.1 or higher
- USB version 2.0 or higher

Tested Bluetooth[®] USB adapters:

- DIGITUS[®] Bluetooth[®] 4.0
- Tiny USB Adapter
- HAMA[®] Nano v2.1

Note

BlueSoleil drivers and software are not supported by IEM GmbH and may lead to complications with device communication.

5.2 Installing on Windows[®]

Procedure:

- Install HMS CS.
- When using a USB cable, install the USB driver.
- Connect the USB cable and/or Bluetooth[®] USB adapter to the computer.

5.2.1 Installing HMS CS from website

1. Download the HMS CS from the following URL: www.iem.de/hmscs.
2. Unzip the downloaded .ZIP file.
3. Click on the Setup.exe file to start the installation program.
4. Select the desired language.
5. Click Next in the new window to start the installation of the HMS CS. The installation wizard appears.
6. Follow the instructions on the screen.
7. This concludes the software installation process.

Installing the USB driver

Note

When using a USB cable, please install the USB driver before connecting the cable to the computer.

1. Click on USB Cable Driver in the installation menu.
2. Follow the instructions displayed on the screen.

This concludes the software installation process.

5.3 Installing on macOS

Procedure:

- Install HMS CS.

5.3.1 Installing HMS CS from website

1. Download the HMS CS from the following URL: www.iem.de/hmscs.
2. Unzip the downloaded .ZIP file.
3. Open the directory macOS.
4. Move the file HMS.dmg to your program directory and execute it.
5. Select the desired language.
6. Click Next in the new window to start the installation of the HMS CS. The installation wizard appears.
7. Follow the instructions displayed on the screen.
8. This concludes the software installation process.

5.4 Updating HMS CS

Updating an already-installed version of Hypertension Management Software CS to the latest version is no different from installing the software from scratch. Existing configurations will not be changed by an update. With regards to GDT (German communication standard) and network settings, the user does not need to do anything when carrying out an update. However, a data backup is recommended before updating Hypertension Management Software Client Server (HMS CS).

We recommend performing HMS CS software updates through our local trained partners.

6 Toolbar

The toolbar is provided at the top of the working window. It contains buttons (icons) used to access important functions quickly.

Tip

If you use the mouse to hover over a symbol, a brief tooltip will appear.

Symbol	Meaning
	New patient
	Patient list
	Device communication
...	Options
	Delete entry
	Export
	Print

7 Starting and exiting HMS CS

Open program

Double-click on the **IEM** icon  on your desktop.

HMS CS is started. Information on the program loading progress is displayed.

Quit program

Click on the **X** at the top right of the working window.

8 First steps with the sample patient

Once you have successfully installed HMS CS, you can try out key functions using the sample patients *John Doe* or *John Doe Jr. (Junior)*.


In the case of the sample patient *John Doe Jr. (Junior)*, the diagnosis also includes information on the study used, which you can configure under global children's blood pressure limits (see also chapter "Global blood pressure limits").

1. Double-click on the **IEM** icon  on your desktop.

HMS CS is started. Information on the loading progress is displayed.

The working window then appears.

2. Click on the **Patient list** icon .

The **Patient list** window then appears.

3. Click on the line with **John Doe**, then on **Open Patient**.

The **Patient information** tab for John Doe is displayed.

The following tabs are available:

- Patient information
- Blood pressure-PWA


8.1 Patient information

The patient information tab features several sections: address, patient data (patient ID, date of birth, weight, etc.), blood pressure limits, emergency contacts, medical history and medication.

8.2 Blood pressure–PWA

1. Click on the **Blood pressure–PWA** tab to display the measurement data for John Doe.

The left side of the **Blood pressure–PWA** tab contains a list of different categories and the previous measurements.

2. For example, click on one of the existing measurements in the 24-hour ABPM category.

The selected measurement turns green and the associated measurement data table is displayed.

The values highlighted in red are measurements that fall outside of the defined limits.

3. To display other analyzes, click on the desired analysis icons.
4. To print out the selected measurement, click on the print icon in the toolbar.

Analysis toolbar:

Tip

If you use the mouse to hover over a symbol, a brief tooltip will appear.

9 Editing patient information

Patient information is stored in a database. You can:

- enter new patients,
- edit data of patients already stored,
- import existing patient information that you have already access to from a different source (GDT).

All patient information can be edited at any time once entered.

9.1 Signing-in a new patient

1. Click on the **New Patient** icon  in the toolbar.

The **New Patient** window appears.

Note

Patient ID, Surname and **Date of birth** are mandatory fields (this information is used for sorting and searching); all other information is optional.

2. To save the new patient, click on **Save**.

To discard the new patient, click on **Cancel**.

Either option will return you back to the working window.

The **Patient information** tab shows the data for the new patient. It contains several sections: address, patient data, blood pressure limits, emergency contacts, medical history and medication.

9.2 Selecting an existing patient

Select a patient from the patients already input into HMS CS to

- edit their patient information,
- view their previous measurements,
- prepare the monitor for this patient or
- import the values from a measurement that has just been performed from the monitor into HMS CS.

1. Click on the **Patient list** icon  in the toolbar.

The **Patient list** window appears with all patients already entered into HMS CS.

2. When you see the desired patient in the window, click on the corresponding list entry and then on **Open Patient**.

Tip

Double-clicking on the list entry using the left mouse button is an even faster way of achieving this.

To search for a patient:

1. Enter the **Surname, Forename** or **Patient ID** at the top right in the search field.

While you are entering the value, HMS CS will dynamically search through the patient list and display the patients found.

2. Click to select the appropriate list entry and then click on **Open Patient**.

If you cannot find the patient you are looking for:

Click on **New Patient** to enter a new patient.

9.3 Adding to and changing patient information

To change an address or patient data, double-click on the field to be changed and enter the new information in the relevant fields.

To save the changes, click on **Save**.

To add **Emergency contacts**, **Medical history** and Medication, click on **New entry** in the relevant field.

Another window entitled **Emergency contact** appears.

1. Enter the new information into the relevant fields.

2. To accept the new information, click on **Save**.

The window will close.

9.4 Defining individual blood pressure limits

Click on the button **ABPM**, **Individual** or **Dipping** in the **Blood pressure limits** field in the patient information tab.

Double-click on the relevant field to open an editing window where you can set the blood pressure limit values for the patient currently open or reset it to the default values.

Change the values as desired, then click on **Save**.

If the limit values are exceeded, the measurement results will be marked accordingly in the analysis.

9.5 Deleting a patient

Open a patient and click on the **Delete** icon  in the toolbar.

Click on **Yes** to confirm that you wish to delete the patient.

The current patient will be removed from the database, including all measurement data.

9.6 Importing blood pressure readings manually

The special functional field **Manual BP** is provided for blood pressure measurements that you as the doctor or the patient with their own blood pressure measurement has taken.

To store your blood pressure measurements in HMS CS, you can

- enter values provided to you by the patient yourself
- directly import values recorded by the blood pressure monitor from the monitor via data transmission

To enter measurements manually:

1. In the working window, click on the **Blood pressure-PWA** tab.
2. Click on **Manual BP** in the list on the left.
The field **Manual BP** is now green.
3. Click on **Manual BP** again, but this time using the right mouse button
4. Select the field **Manual measurement**.
5. Enter the measured values, date and time (current date and time are automatically entered).
6. Click on **Save**.

10 Program configurations (options)

Click on the **Options** icon  in the toolbar.

The following tabs are available:

- General
- Port Settings
- Report
- GDT settings
- Export
- Audit trail
- Customization

10.1 General

The following HMS CS settings can be **viewed** and possibly modified under **General**.

- General
- Units
- Calibration
- Language
- Database
- Global blood pressure limits
- About
- Device names
- Import
- Data back-up

10.1.1 General

Under **General** you can have the patient list displayed upon program launch, activate Bluetooth® and summarize SBPM measurement series using the **Activation switches**.

You can also modify the pause time between individual measurements for triple PWA.

10.1.2 Units

It is possible to change the **Weight unit**, **Length unit** and **Resist unit**.

Click on the field to be changed and select the desired unit.

10.1.3 Calibration

You can use the **Allow external calibration** activation switch to allow external calibration of the monitors and determine the frequency of calibration.

10.1.4 Language

The user interface is available in various languages.

To change the language:

1. Click on the **Options** icon in the toolbar, then on **Language** in the **General** field.

The **Language** window then appears.

2. Select the desired language from the drop-down list.

3. Click on **Save**.

The window will close.

4. To enable the changes to take effect, exit and restart HMS CS.

10.1.5 Database

Patient data and the associated measurement data are stored in a database. This is where settings are made regarding access to the database. Please refer to your specialist supplier or directly to IEM GmbH for further information.

1. Click on the **Options** icon in the toolbar, then on **Database**.

The **Database** window then appears.

2. Make the desired changes.

3. Click on **Save**.

The window will close.

4. To enable the changes to take effect, exit and restart HMS CS.

10.1.6 Global blood pressure limits

You can specify global limit values for systolic and diastolic blood pressure. If the limit values are exceeded, the measurement results will be marked accordingly in the analysis.

These values are automatically saved as limit values with newly signed-in patients.

The blood pressure limits in children and adolescents aged 0 to 17 years are calculated using the reference tables in the KiGGS study¹, which give a comprehensive description of the health situation of children and adolescents living in Germany. This has shown that the blood pressure limits of children and adolescents depend on age, sex and size.

In 2010, the European Society for Hypertension (ESH) published extensive tables², on which the HMS CS limit values are based. The limit values are determined according to the 95% percentile curve.

The limit value is the one that is equal or less for 95% of a collective (statistical analysis of 15,000 children).

All values above that are deemed to be hypertension.

Working with the percentile curve:

To display the percentile curve (only with patients aged between 3 and 17 years), the date of birth must be entered for the patient, as this provides the basis for HMS CS to calculate the patient's age.

The analysis is always based on the patient's current age. To create a patient history, a print-out must be created for each appointment.

Note

- Blood pressure limits can be set individually for each patient under the **Patient information** tab (see also chapter "Defining individual blood pressure limits").

¹ Neuhauser et al.; Blood pressure percentiles by age and height from nonoverweight children and adolescents in Germany. Pediatrics. 2011 Apr;127(4):e978-88. doi: 10.1542/peds.2010-1290. Epub 2011 Mar 7. PMID: 21382947.

² Lurbe et al.; European Society of Hypertension. Management of high blood pressure in children and adolescents: recommendations of the European Society of Hypertension. J Hypertens. 2009 Sep;27(9):1719-42. doi: 10.1097/HJH.0b013e32832f4f6b. PMID: 19625970.

Specifying blood pressure limit values for analysis

1. Click on the Global blood pressure limits button.
The Global blood pressure limits window appears.
2. Enter the desired limit values for adults.

 **Note**

Dipping is the physiological decline in nocturnal blood pressure. The lowest values generally occur between 2 am and 3 am.

Dipping variants:

- Inverted: Nocturnal blood pressure increase in ratio to daytime blood pressure
- Non-Dipper: nocturnal blood pressure decrease in the 1-digit range from 0 to < 10% in ratio to daytime blood pressure
- Normal: Nocturnal blood pressure reduction starting at 10% to < 20% in ratio to daytime blood pressure
- Extreme Dipper: Nocturnal blood pressure reduction of 20 % or more in ratio to daytime blood pressure

3. To accept the new limit values, click on **Save**.

The window will close.

Specifying children's blood pressure limit values for analysis

1. Click on the Global blood pressure limits button.
The Global blood pressure limits window appears.
2. Click on the Children button.
There are various studies to choose from for children.
3. Use the Function button to select the desired study.
The corresponding literature reference is displayed.
4. The category and percentile can be selected on the basis of the study; to do so, click on the desired functional field.

10.1.7 About

Clicking on the **About** Button will display manufacturer information.

The following information about HMS CS is displayed:

- Manufacturer information
- Name and version of the software
- Versions of the program libraries
- Version and date of the last build

10.1.8 Import ABDM.mdb

ABDM.mdb is a database file from our previous hypertension management software (ABPM). This file was used to store all patient information and the associated measurement series. The **ABDM.mdb** button can be used to import this ABPM database into the current HMS CS database.

CAUTION

Please note that when importing ABDM.mdb the existing HMS CS database will be overwritten. Should you be considering importing a previously created database, please contact your specialist supplier or IEM's Technical Support directly.

10.1.9 Import GDT

The GDT button can be used to import practice data and patient information through a corresponding GDT file. Select the GDT file to be imported from the directory and click on **Open**. The GDT file will then be imported; once this is done, click on **Save**.

Note

More information is available in the chapter "GDT settings".

10.1.10 Import patients

The **Patient** button can be used to import previously exported patient data. Select the patient file to be imported from the directory and click on **Open**. The patient file will then be imported; once this is done, click on **Save**.

10.1.11 Data backup

To back up your data, proceed as follows:

1. Click on the **Back up data** button.
The **Back up data** window appears.
2. Enter a filename and location to save the data you are backing up.
3. Click on **Save**.

Your data will be backed up. Depending on the size of the database, this process may take several minutes.

10.1.12 Restore data

To restore your data, proceed as follows:

1. Click on the **Restore data** button.
The **Restore data** window appears.
2. Select the appropriate file containing your backed-up data.
3. Click on **Open**.
4. Confirm that you wish to restore your data.

Your data will be restored. Depending on the size of the database, this process may take several minutes.

 **Note**

Note: When restoring the data, the current database will be overwritten!

10.2 Port settings

This is where you define the interface that will be used to connect the device to the computer. The connection options are:

- Bluetooth[®]
- Cable with serial interface (COM1, COM2...)
- Cable with USB interface

 **Note**

The connection interface to be used is dependent on the device.

10.2.1 Bluetooth® interface

Defining the Bluetooth® interface for the device

1. Click on the **Bluetooth®** tab, then on the button **Add device**.
The **Bluetooth®** message window appears.
2. Switch the monitor on and switch to the device's pairing mode (read the information in the relevant chapter of the individual devices).
3. Click on **OK** in the **Bluetooth®** message window.
After a moment, the serial number of the blood pressure monitor will appear in the window, for instance CP3327.
4. Click on the corresponding **Serial number**, then on the button **Pairing**.
The **Pairing** message window appears, marking the end of the pairing process.

 **Note**

When using the operating system Windows®, the message Add device will appear in the taskbar. Open the window and click on the button Allow. The pairing code for all IEM devices is: 6624.

5. Click on **OK**.
The new device is displayed in the Bluetooth® device list.
6. Click on **Save** in the **Interfaces** window.

10.2.2 Serial/USB interface

Defining a serial/USB interface for a device

1. Click on serial USB, then on Add device.
2. To search for a device, first establish a connection with the PC, then switch on the device (read the information in the relevant chapter of the individual devices).
3. Click on the Search button.
The device that is found will be displayed in the Device connection window (the Interface and Type fields will be automatically completed). If no device is found, a message to this effect will appear (the fields Interface and Type remain empty in this case).
4. Click on Save.
The new device is displayed in the interface list.

Perform connection test for serial/USB interface:

1. Select the device to be tested in the device list.
2. Click on the **Connection test** button.

A corresponding message appears stating whether the connection test was successful.

10.2.3 Deleting the monitor from the list

1. Click on the **device to be removed** in the **Interface** window (serial/USB or Bluetooth).
2. Confirm the query **Really delete?** by clicking on **Yes**.

The device will be removed from the list.

10.3 Report

There are various options available under **Report**.

- Initial: Additional start period in addition to the existing four time periods.
- Static start: Start time of graphic displays
- Report selection: It is possible to select here which data is included in the findings report.
- Axis range: The axis range for the analysis can be set here.

Under **Details** you may choose from the following standards:

- Absolute or relative vascular age
- With or without smiley
- With or without MAP-C2 calibration (PWA)
- Min. and max. time stamps

Activate or deactivate the relevant activation switches.

Once you have made the desired settings, save the changes by clicking on **Save**.

10.4 GDT settings

GDT (GeräteDatenTransfer, device data transfer) is a data exchange format used by private practices in the German healthcare system. The GDT interface is used to enable system-diagnostic data transmission between medical measurement devices and the PMS (**Patient Management Software**).

GDT must be configured to enable the automatic exchange of patient data between your PMS and HMS CS. If the settings are correct, HMS CS can be started from your office software and patient data can be imported directly.

1. Click on the tab **GDT settings**.
2. Under **Settings**, click on the button **Selection**.
3. Here you can specify the directory shared between HMS CS and your PMS. HMS CS and the PMS must be set to this same directory. Ideally, the HMS CS program directory should be set.
4. In the field **PMS -> HMS CS file**, enter the name of the GDT file used to transmit the patient data of your PMS to HMS CS. The same name must be set in HMS CS and in your PMS.
5. In the field **HMS CS -> PMS file**, enter the name of the GDT file used to transmit the findings report from HMS CS to your PMS. The same name must be set in HMS CS and in your PMS.
6. Click on **Save**.
7. Enter the start file **HMS_GDT.exe** into the settings of your PMS.

Exporting Excel, XML and PDF files

1. Click on **GDT settings**.

Under **Export GDT** there is the option

- to select one of the following encoding types:
 - ISO-8859-1
 - IBM437
 - ASCII
- to export currently selected measurement series in the following file formats:
 - XML
 - XLS
 - PDF

2. Click on the **Selection** button for the desired file type and set the storage location and name for the file.

10.5 Export

10.5.1 Export file name

The name of the exported file can be set under the **Export** tab.

1. Click on the first functional field under **Export file name**.
2. Select what the file name should begin with.
3. Proceed with the other functional fields in the same way.
4. To save the changes, click on **Save**.

10.5.2 Formatting

The date and time format can be set here. To do so, click on the desired field and change the value. Several formats are available.

10.5.3 Exporting as CSV

The activation switch allows the following types to be activated for **CSV exports**:

- CSV Export ABDM
- CSV Export PWA

10.6 Audit-Trail

The audit trail is a quality assurance tool that is used to monitor and record changes to processes.

The switch **Audit trail** can be activated so that changes are recorded. You can also search for individual patients and search between individual patient groups and patient information.

10.7 Customization

The following options are available under Customization:

- Print
- Logos
- Colors

10.7.1 Printing

Under **Print** you can modify the **Header** and **Footer** and add or modify an **Office stamp** and/or **Logo for your practice**.

1. To do this, click on the relevant field and modify it as desired.
2. To save the changes, click on **Save**.

10.7.2 Logos

The logo for the application and the background image for HMS CS can be modified here.

1. Click on the **Select** button and select the desired file.
2. Click on **Open**.

The application logo and/or background image will be uploaded.

10.7.3 Colors

The colors for the analyzes (blood pressure measurements) can be defined here.

1. To do so, click on the desired color or line in the relevant functional field.
2. To save the changes, click on Save.

11 Printing

The **print** function allows you to selectively print out analyzes.

1. Click on the **Patient list** icon  in the toolbar and select the required patient and measurement data.
2. To print, click on the **Print** icon  in the toolbar.
3. The **Printing** window then appears.
4. Select which printer to use and which reports to print.
5. The **Save** button allows you to save your settings here.
6. Set the desired page format.
7. Click on the **Print** button.

11.1 Printing in-office blood pressure measurement

For **in-office blood pressure measurements without PWA**, the following printing options are available:

- Print range
- Reports: patient data sheet, findings, profile, bar chart and measurement data table
- Page format
- Save as PDF

For **in-office blood pressure measurements (PWA)**, the following printing options are available:

- Reports: patient data sheet, findings, profile, profile (PWA), bar chart, measurement data table and patient report
- Page format
- Save as PDF

11.2 Printing 24-hour blood pressure measurement

For **24-hour blood pressure measurements (ABPM)**, the following printing options are available:

- Print range
- Reports: patient data sheet, findings, profile, bar chart, measurement data table
- Hourly mean values
- Page format
- Save as PDF

For **24-hour ABPM PWA**, the following printing options are available:

- Reports: patient data sheet, findings, profile, profile (PWA), bar chart, measurement data table, patient report
- Hourly mean values
- Page format
- Save as PDF

12 Exporting measurement data

You can export a patient's complete data or individual measurements.

12.1 Exporting a patient's complete data

1. Click on the **Patient list** icon  in the toolbar and select the required patient. Then click on **Open patient**.

The patient information will be displayed.

2. Click on the **Export** icon .

3. Click on the **Export** button.

A new **Upload** window opens.

4. Select the save location and enter a file name, then click on **Save**.

5. The "Success!" window appears – confirm this by clicking on **OK**.

6. Click on **Close**.

12.2 Exporting individual measurements

1. Click on the **Patient list** icon  in the toolbar and select the required patient and measurement data under **Blood pressure-PWA**.

2. Click on the **Export** icon  in the toolbar.

3. Under **Measurement**, select the file type and click on the **Export** button.

The **Export measurement series** window opens.

4. Select the save location and enter a file name, then click on **Save**.

5. The "Success!" window appears – confirm this by clicking on **OK**.

6. Click on **Close**.

13 Data analysis

Tip

If you use the mouse to hover over a symbol, a brief tooltip will appear.

The following analyzes and functions for the analysis of measurements are available under the tabs **Blood pressure–PWA**:

Blood Pressure – PWA

Measurement series

Icon	Description
	Table of measurements
	Profile
	Bar chart
	Report
	Hourly mean values
	Profile (PWA)
	Patient report

Single measurement

Icon	Description
	Pulse wave analysis
	Arterial stiffness
	Patient report

13.1 Measurement data analysis

Pre-condition:

- The measurements have been uploaded from the blood pressure monitor and are stored in HMS CS

1. Select the required patient.
2. In the working window, click on the Blood pressure–PWA tab.

The left side of the Blood pressure–PWA tab contains a list of previously performed measurements.

3. Click on a measurement date.
The associated measurement data table is displayed.
The values highlighted in red are measurements that fall outside of the defined limits.
4. To display additional analyzes, click on the required analysis icon.

13.2 Entering findings for measurement series

1. Click on the measurement date with the right mouse button. Select the option Findings in the context menu.

The Measurement series window appears.

2. Enter your findings in the field Findings/Comment.
3. To save the findings, click on Save.

The window will close.

13.3 Blood pressure – PWA tab

13.3.1 Table of measurements

The **Measurement data table**  icon lists all measurement values of a measurement series in table format.

To display the measurement data table, click on the **Blood pressure–PWA** tab, then on the **Measurement data table** analysis icon .

The values highlighted in red are measurements that fall outside of the defined limits.

The comments column is used to automatically save special events such as the EVENT button being pressed or an error message.

Commenting on measurement values

1. Click on the **Comment** column in the desired line.
2. Enter your note.
3. Press the enter button.

Excluding measurements

If one measurement is completely inconsistent and implausible and risks distorting a representative long-term analysis, it can be excluded:

1. Click on the right mouse button and deactivate the activation switch.
2. Click on the left mouse button to deactivate the corresponding column.
The measurement number disappears and the measurement is henceforth excluded from statistical analyzes.
3. To re-include the measurement, click on the right mouse button, then on the activation switch.
4. Click on the left mouse button to re-activate the corresponding column.

Printing the measurement data table

Click on the **Print** icon  in the toolbar.

13.3.2 Blood pressure profile

This profile graph shows the following values in the measurement series as a line chart:

- Systolic values
- Diastolic values
- Blood pressure limits
- Heart rate
- Average values.

To display the profile, click on the **Blood pressure-PWA** tab, then on the **Profile** analysis icon .

The left y-axis with the unit mmHg represents the systolic, diastolic and average values (blood pressure values). The right y-axis with the unit bpm represents the heart rate.

The x-axis represents time. The four adjustable time periods are visually highlighted.

The upper blood pressure limits (systolic, diastolic) are displayed as horizontal set point curves.

The night-time interval is shown in shaded grey, with the moon symbol marking the start and the sun symbol marking the end.

Note

Individual blood pressure limits:

These values are defined under the **Patient information** tab in the **Blood pressure limits** section.

Global blood pressure limits:

These values are defined in **Options** under the **General** tab in the **Global blood pressure limits** section.

Displaying and hiding individual functions:

Heart rate

Click on the activation switch **Hr** (heart rate) to display or hide.

Average values

Click on the activation switch **MAP** (mean arterial pressure) to display or hide.

Display and hide hourly mean values

To modify the desired number of hours for the hourly average value, click on the functional field **Hourly mean values** and select the desired average. (The tab may need to be reloaded.)

Battery voltage

Click on the **Battery voltage** activation switch. The battery charge level over time is displayed as a 24-hour curve parallel to the blood pressure.

Individual values

Click on the diagram with the left mouse button.

A vertical line appears and the measured values are displayed in a new window.

To deactivate the display again, move the mouse outside of the diagram or click again with the left mouse button.

Tip

To see adjacent measurement values, move the mouse over the diagram. The vertical line follows the movement of the mouse and the respective measurements are displayed.

Magnify diagram range

Click on the diagram with the left mouse button and hold the mouse button down. Now drag the frame down around the area to be magnified from **left to right**, then release the mouse button.

Restoring original diagram size

Click on the diagram with the left mouse button and hold the mouse button down. Now drag a line from **right to left** before releasing the mouse button.

13.3.3 Bar chart

This profile graph shows the following values in the measurement series as a bar chart:

- Systolic values
- Diastolic values
- Blood pressure limits
- Heart rate

To display the bar chart, click on the **Blood pressure–PWA** tab, then on the **Bar chart** analysis icon .

The left y-axis with the unit mmHg represents the systolic, diastolic and average values (blood pressure values). The right y-axis with the unit bpm represents the heart rate.

The x-axis represents time. The time periods are visually highlighted.

The upper blood pressure limits (systolic, diastolic) are displayed as horizontal set point curves.

Note

Individual blood pressure limits:

These values are defined under the **Patient information** tab in the **Blood pressure limits** section.

Global blood pressure limits:

These values are defined in **Options** under the **General** tab in the **Global blood pressure limits** section.

Enlarge diagram range

Click on the diagram with the left mouse button and hold the mouse button down. Now drag the frame down around the area to be magnified from **left to right**, then release the mouse button.

Restoring original diagram size

Click on the diagram with the left mouse button and hold the mouse button down. Now drag a line from **right to left** before releasing the mouse button.

13.3.4 Report

The findings report contains important statistical statements on systolic and diastolic blood pressure. The values are displayed for day, night and total periods, depending on the measurement series selected.

To display the findings report, click on the **Blood pressure–PWA** tab, then on the **Findings report** analysis icon .

Mean values are provided for the measurements directly under day and night.

The mean blood pressure limits are specified here as target values and can be defined under patient information in the blood pressure limits section.

Other target values are specified by HMS CS and cannot be set individually.

Note

Individual blood pressure limits:

These values are defined under the **Patient information** tab in the **Blood pressure limits** section.

Global blood pressure limits:

These values are defined in **Options** under the **General** tab in the **Global blood pressure limits** section.

13.3.5 Hourly mean values

This analysis lists all hourly mean blood pressure and heart rate values in a table.

To display the hourly mean values, click on the **Blood pressure–PWA** tab, then on the **Hourly mean values** analysis icon .

Modifying calculation basis for hourly mean values

Click on the desired number of hours (1, 2, 3, 4, 6, 8) in the functional field **Mean calculation basis (h)**.

The time periods in the left column **Time** will be displayed accordingly. The average hourly value is recalculated.

13.4 Comparing multiple measurement series

If two or more measurement series are stored for a patient, these can be compared against one another.

Depending on the analysis

- diagrams for each of the measurement series will be arranged one above the other or
- the values will be aggregated then displayed as a graph.

Selecting and comparing multiple measurement series

1. Click on the first measurement.

The measurement series will be highlighted.

2. Hold down the **CTRL key**, then click on the other needed measurement series.

These measurement series will also be highlighted.

3. Click on the desired analysis icon.

13.4.1 24h PWA

Note

24-hour PWA analysis is only possible with a Mobil-O-Graph® blood pressure monitor and a license key. Please contact the manufacturer or your specialist supplier with any questions you may have.

This analysis shows you the results of the PWA over a 24-hour period. The following values in the measurement series are also shown alongside the blood pressure readings and heart rate readings in a diagram:

- Central blood pressure (cBP)
- Pulse wave velocity (PWV)
- Cardiac output (CO)
- Peripheral resistance (TVR)
- Heart rate (Hr)

To show the changes in the aforementioned values over time, click on the **Blood pressure – PWA** tab, then on the **Profile (PWA)** analysis icon

13.4.2 Patient report

In the case of PWA measurements as part of 24h ABPM values, the patient report displays the average of selected hemodynamic information determined during the measurement period.

To open the patient report (day average), click on the **Blood pressure–PWA** tab, then on the **Patient report** icon .

The patient report tab is divided into four sections:

Peripheral blood pressure:

The peripheral blood pressure section displays the peripheral systolic and peripheral diastolic blood pressure measurements.

Central blood pressure:

The central blood pressure section displays the calculated central systolic blood pressure.

Vascular stiffness (In the USA: In patients 40 years of age and older):

The vascular stiffness section shows the calculated pulse wave velocity (PWV) and the augmentation index at a heart rate of 75 bpm (AIx@75), which provides an indication of the stiffness of large and small arteries.

Vascular age (Not available in the USA):

The vascular age section shows the vascular age of the patient determined on the basis of hemodynamic information. The vascular age display can be changed in the HMS CS options. The display of the vascular age can be configured to be an absolute or relative value.

Central blood pressure and calibration:

Current technology requires calibration to calculate central aortic systolic blood pressure. Peripheral systolic blood pressure is usually used for this purpose. Calibration can also be performed on the basis of the measured MAP (mean arterial pressure, corresponds to the maximum oscillometric amplitude).

Scientific findings recently published have proven for the first time that calibration based on measured MAP values is more predictive than other comparable methods.³⁴⁵

³ Kazuaki Negishi et al. Importance of Calibration Method in Central Blood Pressure for Cardiac Structural Abnormalities. American Journal of Hypertension 2016

⁴ Athanase D. Protopsaltis et al. Left-ventricular hypertrophy is associated better with 24-h aortic pressure than 24-h brachial pressure in hypertensive patients: the SAFAR study. American Journal of Hypertension 2014

⁵ Thomas Weber et al. Validation of a Brachial Cuff-Based Method for Estimating Central Systolic Blood Pressure. Journal of the American Heart Associations 2011

Further information on the vascular age display:

The vascular age of a patient is calculated based on the pulse wave velocity measured. In addition to the normal age-related development of arterial stiffness in healthy people (green curve), measured in meters per second (m/s), there are also people who do not fit the normal pattern (red dot). If you follow the measurement horizontally along to the normal curve, you can derive the vascular age of a patient.


Fig. 1: PWV by patient age

The vascular age display can be modified in the HMS CS options.

Here, you can define the standard methods for:

- Absolute or relative vascular age
- Patient report display with or without smiley

13.4.3 Arterial stiffness

The arterial stiffness section displays the pulse wave velocity (PWV) (Not available in the USA) measured for large vessels and the augmentation index at a heart rate of 75 bpm (AIx@75) in the case of small vessels.

To show the changes in the aforementioned values over time, click on the **Blood pressure–PWA** tab, then on the **Arterial stiffness** analysis icon .

13.5 ECG tab

The ECG tab is no longer available because the ECG functionality is no longer supported.

The following applies to your old ECG data:

When updating to HMS CS 6.2, a patient's ECG data is still stored in the database. You can therefore continue to access this data, even if it is no longer visible in HMS CS 6.2 or newer:

1. When opening the patient with a previous HMS CS version (HMS CS 6.1 or older), the data is still available and editable.
2. Even in HMS CS 6.2 or newer, you can access the data as follows:
 - Open the desired patient with stored ECG data.
 - Click on the **Export** button and select **ECG&Weight** in the **Patient** area:

- Click on the corresponding **Export** button and select a storage location for the PDF file with the ECG data contained.
- A PDF file will be generated containing all ECG data stored for the patient in the database. The ECG data are presented according to the visualization in the ECG tab from previous HMS CS versions.

13.6 Weight tab

The Weight tab is no longer available because the visualization of Weight data is no longer supported.

The following applies to your old Weight data:

When updating to HMS CS 6.2, a patient's Weight data is still stored in the database. You can therefore continue to access this data, even if it is no longer visible in HMS CS 6.2 or newer:

1. When opening the patient with a previous HMS CS version (HMS CS 6.1 or older), the data is still available and editable.
2. Even in HMS CS 6.2 or newer, you can access the data as follows:
 - Open the desired patient with stored Weight data.
 - Click on the **Export** button and select **ECG&Weight** in the **Patient** area:

- Click on the corresponding **Export** button and select a storage location for the PDF file with the Weight data contained.
- A PDF file will be generated containing all Weight data stored for the patient in the database. The Weight data are presented according to the visualization in the Weight tab from previous HMS CS versions.

14 HMS CS & Mobil-O-Graph®

14.1 Bluetooth® connection

The following procedure applies when using the blood pressure monitor with a Bluetooth® connection.

Note

Bluetooth® is not supported on operating system macOS.

Prepare and perform the long-term measurement

1. Configure the blood pressure monitor and HMS CS software
2. Prepare the blood pressure monitor for long-term measurements
3. Start the 24-hour measurement

Transferring and analyzing the long-term measurement

4. Transferring and storing the long-term measurement results from the blood pressure monitor

14.1.1 Pairing the blood pressure monitor with HMS CS

Pre-conditions:

- Important: computer must support Bluetooth®
- Bluetooth® must be enabled in HMS CS, see chapter "General"
- Computer must be switched on

Note

To verify whether Bluetooth® is enabled in HMS CS, click on the Options icon  in the toolbar and check under General.

Procedure on computer:

1. Start **HMS CS**.
2. Click on the **Patient list** icon  in the toolbar and select the required patient in the **Patient list** window.
3. Click on the **Options** icon  in the toolbar, then on the **Interfaces** tab.
4. Click on the **Bluetooth®** tab.
5. Click on **Add device**.

A window with the following message appears:

"Switch on the device and change to pairing mode."

Then click on 'OK'. Details on accessing pairing mode are provided in the manual."

Procedure on blood pressure monitor:

6. Switch on blood pressure monitor

Switch to Pairing mode:

7. Hold down the START button  and press the DAY/NIGHT button  once, then release both buttons. **bt** will appear flashing on the display.
8. Keep pressing the START button  until the letters **PAlr** flash on the display.
9. Press the EVENT button .

The letters **PAlr** cease flashing and an audible signal is emitted.

Procedure on computer:

10. Click on **OK**.

The **Bluetooth® device** search window appears.

After a moment the serial number appears in the window, for instance C00607.

11. Click on the serial number.

12. Click on **Pairing**.

**Note**

When using the operating system Windows®, the message add device will appear in the taskbar. Click on the message, then click on the button Allow. The pairing code for all IEM devices is: 6624.

The following message appears in HMS CS:

"The (one-time) pairing process was successful"

13. Click on **OK**.

14. Click on **Save**.

The Bluetooth® interface between the blood pressure monitor and HMS CS is now configured and HMS CS will henceforth detect the blood pressure monitor as soon as you switch it to "**bt**" connection mode.

14.1.2 Prepare the blood pressure monitor for long-term measurements

Pre-conditions:

- The blood pressure monitor is switched off
- The computer is switched on
- The interface for the blood pressure monitor is known by HMS CS
- HMS CS is started
- The Bluetooth® connection is active

Note

To verify whether Bluetooth® is enabled in HMS CS, click on the **Options** icon  in the toolbar and check under **General**.

Procedure on computer:

1. Click on the **Patient list** icon  in the toolbar and select the required patient in the **Patient list** window.

Procedure on blood pressure monitor:

2. Switch on the blood pressure monitor.
3. Hold down the EVENT button  for around 4 seconds.

The letters **bt** will appear on the display and an audible signal will be emitted.

Procedure on computer:

The **Device action** window will appear on the computer with the buttons **Prepare device**, **Read-out of values**, **PWA**, **Measurement Triple PWA Measurement** and **Cancel**.

4. Click on **Prepare device**.

The **Prepare device** window appears.

CAUTION

If the battery charge level of the monitor is not adequate for long-term measurement, a message to this effect will be displayed. Please heed this information, as having too low a battery charge level may result in the device ceasing to work!

5. Adjust the log parameters as desired, see also chapter "Configuring the log".
6. If the buttons in the **Prepare device** window are shown in red, click on them.
The corresponding settings are adjusted.
7. Save the settings by clicking on **Save**.
8. Click on **OK** in the confirmation window.
9. HMS CS can be closed now.

**Note**

The blood pressure monitor buzzer sounds and the letters **bt end** will appear first on the display, followed by the time.

14.1.3 Start the 24-hour measurement

Procedure on blood pressure monitor:

1. Set up the monitor on the patient (apply the cuff and connect it to the monitor).

**WARNING**

Please note the warnings in the instructions for use of the blood pressure monitor.

2. Switch on the blood pressure monitor.

The time as set appears on the display.

3. Press the START button  to perform a manual measurement to ensure that the monitor is working as desired.

**Note**

A successful measurement is required to activate the log!

4. If everything is OK, the patient can be released.

14.1.4 Transferring and storing long-term measurement results

Pre-conditions:

- Blood pressure monitor is switched on
- The computer is switched on
- Interface for blood pressure monitor is known by HMS CS

Procedure on computer:

1. Start **HMS CS**.
2. Check whether **Bluetooth[®]** is activated in HMS CS (in the **Options** under **General**).

Procedure on blood pressure monitor:

3. Remove the monitor from the patient (detach the cuff and disconnect it from the monitor).
4. Hold down the EVENT button  for around 4 seconds.

The letters **bt** will appear on the display and an audible signal will be emitted.

Procedure on computer:

The **Device action** window will appear on the computer with the buttons **Prepare device, Read-out of values, PWA, Triple PWA** and **Cancel**.

5. Click on the **Read-out of values** button.

The Read-out of values window appears:

"The Patient ID set in the device (9999999999) refers to John Doe (02/08/45). Should the measurement series be assigned to this patient?"

6. Click on **Yes** if the patient you want is displayed.

The **Measurement series** window appears.

7. If necessary, change the times for the daytime and night-time interval.

8. Document your findings as necessary and amend the other information.

9. Then click on **Save**.

The **Delete measurements** window appears:

"Should The Patient ID and the measurements be deleted from the device?"

ATTENTION: If you do not delete the measurements, no new measurement series will be started, but the next measurements will be added to the existing series"

10. Click on Yes to delete the measurement results from the device or click on **No** to keep the measurement results on the blood pressure monitor.

Data transmission is completed.

 **Note**

Normally the measurement results on the blood pressure monitor are deleted once the results have been transferred. If the blood pressure monitor is prepared for a "new" patient, HMS CS will notify the user if any measurement results from a previous patient are still on the blood pressure monitor.

Procedure on blood pressure monitor:

11. Switch off the blood pressure monitor.

14.2 Connection via cable

Cable connection options:

- via serial interface (COM1, COM2, etc.) or
- via USB interface.

The following procedure applies when using the blood pressure monitor with a cable connection:

Prepare and perform the long-term measurement

1. Connect the blood pressure monitor to the computer
2. Configure the blood pressure monitor in HMS CS
3. Prepare the blood pressure monitor for long-term measurements
4. Start the 24-hour measurement

Transferring and analyzing the long-term measurement

5. Connect the blood pressure monitor to the computer again
6. Transferring and storing the long-term measurement results from the blood pressure monitor


14.2.1 Connecting the blood pressure monitor to the computer via cable

Procedure on Computer:

Cable with serial interface (COM1, COM2...)	Cable with USB interface
1. Connect the cable to the serial interface (COM1, COM2, etc.)	1. Connect the cable to the USB socket.

Procedure on blood pressure monitor:

2. Switch off the blood pressure monitor.
3. Insert the plug into the data socket on the left side of the casing until it clicks into place.

CAUTION

The red point on the plug must align with the red point on the data socket. Do not apply force!

4. Switch on the blood pressure monitor

The letters **co** appear on the display.

14.2.2 Configure the blood pressure monitor in HMS CS

Pre-conditions:

- Blood pressure monitor is switched on
- The computer is switched on
- The blood pressure monitor is connected to the computer.

Procedure on computer:

1. Start **HMS CS**.
2. Click on the **Patient list** icon  in the toolbar and select the required patient in the **Patient list** window.
3. Click on the **Options** icon  in the toolbar, then on the **Interfaces** tab.
4. Click on the **Serial/IR/USB** tab.
5. To search for a device, click on **Add device**.
The **Device connection** window appears.
6. Click on **Search**.
The device found is displayed in the **Type** functional field, while the corresponding interface is displayed in the **Interface** functional field. If no device is found, a message to this effect will appear.
7. Click on **Save**.
8. The new device is displayed in the interface list.
9. Perform a connection test.

14.2.3 Prepare the blood pressure monitor for long-term measurements

Pre-conditions:

- The blood pressure monitor is connected to the computer
- The blood pressure monitor and computer are switched on
- The interface for the blood pressure monitor is known by HMS CS
- The patient is selected

Note

Always use fully charged batteries for a new measurement. Ensure the correct polarity when inserting the batteries.

Procedure on computer:

1. Click on the **Device communication** icon  in the toolbar.
2. Click on **Prepare device** in the new window.

The **Prepare device** window appears.

**CAUTION**

If the battery charge level of the monitor is not adequate for long-term measurement, a message to this effect will be displayed. Please heed this information, as having too low a battery charge level may result in the device ceasing to work!

3. Adjust the log parameters as desired, see also chapter "Configuring the log".
4. If the buttons in the **Prepare device** window are shown in red, click on them.
The corresponding settings are adjusted.
5. Save the settings by clicking on **Save**.
6. Click on **OK** in the confirmation window.
7. You can now close HMS CS.

Procedure on blood pressure monitor:

8. Switch off the blood pressure monitor.
9. Disconnect the cable (pull the plug out of the data socket).

14.2.4 Start the 24-hour measurement

Pre-conditions:

- The blood pressure monitor has been disconnected from the computer.

Procedure on blood pressure monitor:

1. Set up the monitor on the patient (apply the cuff and connect it to the monitor).

**WARNING**

Please note the warnings in the instructions for use of the blood pressure monitor.

2. Switch on the blood pressure monitor.
3. The time as set appears on the display.
4. Press the START button  to perform a manual measurement to ensure that the monitor is working as desired.

**Note**

A successful measurement is required to activate the log!

5. If everything is OK, the patient can be released.

14.2.5 Connect the blood pressure monitor to the computer again

After long-term measurement, transfer the data from the monitor to HMS CS.

Pre-conditions:

- The blood pressure monitor has been fitted to the patient's arm and is switched on

Procedure on blood pressure monitor:

1. Switch off the blood pressure monitor.
2. Remove the monitor from the patient (detach the cuff and disconnect it from the monitor).

WARNING

Please note the warnings in the instructions for use of the blood pressure monitor.

Procedure on computer:

Cable with serial interface (COM1, COM2...)	Cable with USB interface
3. Connect the cable to the serial interface (COM1, COM2, etc.)	3. Connect the cable to the USB socket.

Procedure on blood pressure monitor:

4. Insert the plug into the data socket on the left side of the casing until it clicks into place.

CAUTION

The red point on the plug must align with the red point on the data socket. Do not apply force!

5. Switch on the blood pressure monitor.

The letters **co** appear on the display.

14.2.6 Transferring and storing long-term measurement results

Pre-conditions:

- The blood pressure monitor and computer are switched on
- The blood pressure monitor is connected to the computer
- Interface for blood pressure monitor is known by HMS CS

Procedure on computer:

1. Start **HMS CS**.
2. Click on the **Patient list** icon  in the toolbar and select the required patient or sign-in a new patient.
3. Click on the **Device communication** icon  in the toolbar.
4. Click on **Read-out device** in the new window.

The **Read-out of values** window appears:

"The Patient ID set in the device (9999999999) refers to John Doe (02/08/45). Should the measurement series be assigned to this patient?"

Should this series of measurements be attributed to this patient?"

5. Click on **Yes** if the patient you want is displayed.

The **Measurement series** window appears.

6. If necessary, change the times for the daytime and night-time interval.
7. Document your findings as necessary and amend the other information. Then click on **Save**.

The **Delete measurements** window appears:

"Should The Patient ID and the measurements be deleted from the device?

ATTENTION: If you do not delete the measurements, no new measurement series will be started, but the next measurements will be added to the existing series."

8. Click on **Yes** to delete the measurement results or click on **No** to keep the measurement results on the blood pressure monitor.

Data transmission is completed.

 **Note**

Normally the measurement results on the blood pressure monitor are deleted once the results have been transferred. If the blood pressure monitor is prepared for a "new" patient, HMS CS will notify the user if any measurement results from a previous patient are still on the blood pressure monitor.

Procedure on blood pressure monitor:

9. Switch off the blood pressure monitor.
10. Disconnect the cable (pull the plug out of the data socket).

15 HMS CS & Tel-O-Graph[®] BT

15.1 Connection via Bluetooth[®]

The following procedure applies when using the blood pressure measuring device with a Bluetooth[®] connection:

 **Note**

Bluetooth[®] is not supported on operating system macOS.

Prepare and perform the blood pressure measurement

1. Pair the blood pressure measuring device with HMS CS
2. Prepare the blood pressure measuring device for blood pressure measurement
3. Start blood pressure measurement

Transfer and analysis of blood pressure measurement

4. Transfer and storage of blood pressure measurement results from the blood pressure measuring device

15.1.1 Pairing the blood pressure measuring device with HMS CS

Pre-conditions:

- Bluetooth[®] must be enabled in HMS CS, see chapter "General"
- The computer is switched on

 **Note**

To verify whether Bluetooth[®] is enabled in HMS CS, click on the **Options** icon  in the toolbar and check under **General**.

Procedure on computer:

1. Start **HMS CS**.
2. Click on the **Patient list** icon  in the toolbar and select the required patient in the **Patient list** window.
3. Click on the **Options** icon  in the toolbar, then on the **Port settings** tab.
4. In the **Port settings** window, click on the **Bluetooth[®]** tab.
5. Click on **Add device**.

A window with the following message appears:

"Switch on the device and change to pairing mode.

Then click on 'OK'. Details on accessing pairing mode are provided in the manual."

Procedure on the blood pressure measuring device:

6. Switch the blood pressure measuring device on using the  button.

Switch to **Pairing mode**:

7. Hold down the  button until you hear a second audible signal, then release the button.
Wait until **PAI P** flashes on the display.

8. Press the  button again.

An audible signal sounds and **PAI P** appears on the display, this time not flashing.

**Note**

After 3 seconds, the display will show **bt**; ignore this and hold down the button for another 3 seconds.

After 6 seconds, the menu automatically opens and the display will automatically show the next menu item. The display order is:

- **Passive pairing (PAI P)**
- Infrared transmission (ir)
- Active pairing (PAI A)
- Bluetooth[®] transmission (bt)
- Delete measurements (c lr)

Procedure on computer:

9. Click on **OK**. The **Searching for Bluetooth[®] devices** search window appears.

After a moment the serial number appears in the window, for instance T80003T2.

10. Click on the serial number.

11. Click on **Pairing**.

**Note**

When using the operating system Windows[®], the message **Add device** will appear in the taskbar. Click on the message, then click on the button **Allow**. The pairing code for all IEM devices is: 6624.

The following message appears:

"The (one-time) pairing process was successful"

12. Click on **OK**.

13. Click on **Save**.

The Bluetooth[®] interface between the blood pressure measuring device and HMS CS is now configured and HMS CS will henceforth detect the blood pressure measuring device as soon as you switch it to "**bt**" connection mode.

15.1.2 Prepare the blood pressure measuring device for blood pressure measurement

Pre-conditions:

- The blood pressure measuring device is switched off
- The computer is switched on
- The interface for the blood pressure monitor is known by HMS CS
- HMS CS is started
- The Bluetooth[®] connection is active

 **Note**

To verify whether Bluetooth[®] is enabled, click on the **Options** icon  in the toolbar and check under general.

Procedure on computer:

1. Click on the **Patient list** icon  in the toolbar and select the required patient in the **Patient list** window.

Procedure on the blood pressure measuring device:

2. Switch the blood pressure measuring device on using the  button.
3. Hold the  button for 3 seconds until the **bt** appears flashing on the display.
An audible signal sounds and **bt** appears on the display, this time not flashing.

Procedure on computer:

The **Device action** window will appear on the computer with the buttons **Prepare device**, **Read-out of values**, **PWA Measurement**, **Triple PWA Measurement** and **Cancel**.

4. Click on **Prepare device**.
The **Prepare device** window appears.
5. Modify the configuration as you desire.
6. If the buttons in the **Prepare device** window are shown in red, click on them.
The corresponding settings are adjusted.
7. Save the settings by clicking on **Save**.
8. You can now close HMS CS.

 **Note**

You will hear the blood pressure measuring device buzz and the start screen will appear on the display.

15.1.3 Start blood pressure measurement

Pre-conditions:

- The blood pressure measuring device is switched on

Procedure on the blood pressure measuring device:

1. Set up the monitor on the patient (apply the cuff and connect it to the monitor).

**WARNING**

Please note the warnings in the instructions for use of the blood pressure measuring device.

2. Press the button to start the measurement.

The blood pressure measuring device confirms this with a brief audible signal and briefly shows the display functions. The cuff slowly pumps up. The current pressure is shown on the display. Once the initial pumping process is complete, a second pumping process will occur. Once the monitor detects a pulse, the icon appears. The blood pressure measuring device will pump up the cuff once again for PWA measurement if PWA has been activated in HMS CS. The measuring process occurs as the air is let out. The blood pressure measuring device confirms that the measurement is complete by emitting a brief audible signal.

3. Wait for the measurement.

**Note**

A measurement must have been completed successfully before it can be transmitted to HMS CS.

15.1.4 Transferring and storing blood pressure measurement results

Pre-conditions:

- The blood pressure measuring device is switched on
- The computer is switched on
- The interface for the blood pressure measuring device is known by HMS CS

Procedure on computer:

1. Start **HMS CS**.
2. Check whether Bluetooth® is activated in HMS CS (in the **Options** under **General**).

Procedure on the blood pressure measuring device:

3. Remove the monitor from the patient (detach the cuff and disconnect it from the monitor).
4. Hold the  button for 3 seconds until the **bt** appears flashing on the display.
An audible signal sounds and **bt** appears on the display, this time not flashing.

Procedure on computer:

The **Device action** window will appear on the computer with the buttons **Prepare device**, **Read-out of values**, **PWA Measurement**, **Triple PWA Measurement** and **Cancel**.

5. Click on **Read-out of values**.

The **Read-out of values** window appears.

"The Patient ID set in the device (9999999999) refers to John Doe (02/08/45). Should the measurement series be assigned to this patient?"

6. Click on Yes if the patient you want is displayed.

The **Measurement series** window appears.

7. Document your findings as necessary and amend the other information.

8. Click on **Save**.

The Delete measurements window appears:

"Should The Patient ID and the measurements be deleted from the device?"

ATTENTION: If you do not delete the measurements, no new measurement series will be started, but the next measurements will be added to the existing series."

9. Click on **Yes** to delete the measurement results or click on **No** to keep the measurement results on the blood pressure monitor.

Data transmission is completed.

**Note**

Normally the measurement results on the blood pressure measuring device are deleted once the results have been transferred. If the blood pressure measuring device is prepared for a "new" patient, HMS CS will notify the user if any measurement results from a previous patient are still on the blood pressure monitor.

Procedure on the blood pressure measuring device:

The blood pressure measuring device will switch off automatically.

16 Prepare device

Pre-conditions:

- Batteries are inserted
- The device is switched on
- The computer is switched on and HMS CS has been opened
- The device is configured in HMS CS (paired)

1. Select the required patient.
2. Connect the device with HMS CS via Bluetooth® (cable).
The **Device action** window opens.
3. Once the connection with HMS CS has been established, click on the **Prepare device** button.
The **Prepare device** window opens.

16.1 Patient list

You can use the patient list to select another patient.

1. Click on patient list in the Prepare monitor window.
2. Select the required patient or sign-in a new patient.

16.2 Setting device clock

The monitor's internal clock will be synchronized with that of the computer.

1. Click on set device clock in the Prepare monitor window.
2. Confirm by clicking on **OK**.

The synchronized time appears on the monitor.

16.3 Importing patient ID

The patient ID of the selected patient is saved in the monitor. Later, when uploading the long-term measurements, this will allow HMS CS to recognize the patient automatically.

1. Click on Send patient ID in the **Prepare device** window.
2. Confirm by clicking on **OK**.

16.4 Testing the device

Ensure that the monitor is working correctly.

1. Click on Device tests... in the **Prepare device** window.

The **Device tests** window appears.

2. Click on the appropriate buttons.

The following tests are available:

Mobil-O-Graph®

Display, keypad, buzzer, version, memo voltage, battery, serial number, calibration date and PWA status

Tel-O-Graph®

Buzzer, version, memo voltage, keypad, battery, serial number, calibration date and PWA status

3. If asked to confirm, click on **OK**.
4. End the test by clicking on **Close**.

16.5 Deleting old measurements/records

Measurements/records are normally deleted in the device when uploaded to the computer. If "old" measurements/records from a previous patient are still present in the monitor, the HMS CS software will advise you of this when preparing the device for a "new" patient.

To remove the "old" measurements/recordings from the device, click on **Delete measurements/recordings** in the **Prepare device** window.

Confirm by clicking on **Yes**.

16.6 Special Mobil-O-Graph® settings

16.6.1 Configuring the log

Click on the desired day period in the Monitor window.

Under day period, set the following:

- the times (start of time period),
- the number of measurements within the period,
- whether the measurements are displayed on the blood pressure monitor (measurement display),
- whether an audible signal (buzzer) is emitted during measurement and
- the number of PWA measurements within the period.

Note

24-hour PWA measurement is only possible with a license key, Mobil-O-Graph® and Bluetooth® interface; please direct any queries to the manufacturer.

16.6.2 Pre-set

The pre-set function lets you save the measurement log that you want.

1. Configure the log that you want in the **Prepare device** window.
2. Click on the **Pre-set** list option.
3. Enter the desired configuration name.
4. Click on **Save**.

When next preparing the monitor, the pre-set measurement log will be available to use. Click on the pre-set list option to open the saved measurement log.

16.6.3 In-office blood pressure monitoring

In-office blood pressure monitoring aims to reduce the workload of the office personnel, improve the quality of treatment and enhance patient comfort. The blood pressure monitor can be carried by the patient around the office (e.g. in the waiting room) and the measurement series is transmitted directly to an office computer via Bluetooth®. Each measurement is sent directly, wirelessly and automatically to HMS CS where it can be analyzed directly by the doctor.

In-office monitoring can be used:

- to create a detailed & accurate summary profile of the patient

WARNING

This system must not be used to provide blood pressure monitoring services with alarm signals in operating theatres or intensive care units.

Preparation of blood pressure monitor for in-office blood pressure monitoring

For in-office monitoring the Bluetooth® interface of the blood pressure monitor is used. If you have not used the Bluetooth® interface before now, follow the instructions under the chapter "Bluetooth® connection".

Pre-conditions:

- The blood pressure monitor and computer are switched on
- The interface for the blood pressure monitor is known by HMS CS (paired)
- The blood pressure monitor is connected with HMS CS via Bluetooth®

1. Once the Bluetooth® connection between HMS CS and the blood pressure monitor has been established, click on the **Prepare device** button.

The **Prepare device** window opens.

2. Activate the **In-office** and **Bluetooth®** switches.
3. Also activate the PWA switch if necessary; this requires a PWA license.
4. Set the desired time periods. Options here include 30,15,12, etc.
5. If the buttons in the prepare monitor window are shown in red, click on them.
The corresponding settings are adjusted.
6. Save the settings by clicking on **Save**.
7. Position the cuff on the patient and connect the hose of the cuff to the blood pressure monitor.

WARNING

Please note the warnings in the instructions for use of the blood pressure monitor.

8. Ensure that the blood pressure monitor is working as desired, then press the START button  to start a manual measurement.
9. Wait for the first automatic measurement to be completed, then check whether the measurement results have been transmitted to HMS CS.

Note

A successful measurement is required to activate in-office measurements!

Assigning the received measurement series

After the first measurement, the **In-office monitoring** window will appear in HMS CS.

10. Click on **Assign**.

The **Selection** window appears.

"Please select a patient to assign the measurement series to."

11. The measurement series can be assigned to either the patient currently open or to another patient in the patient list.

17 Pulse wave analysis

HMS CS in conjunction with the blood pressure monitor offers not only conventional 24-hour blood pressure measurement but also an integrated pulse wave analysis (PWA) system for in-office measurement. This function can be enabled with a PWA license key or a PWA license dongle. The license key / dongle can be obtained from IEM GmbH or your specialist supplier.

Pulse wave analysis is based on the concept that the arterial blood pressure curve contains hemodynamic information that goes beyond peripherally measured blood pressure values. This is used to analyze all information related to the central aortal pulse wave.

The following values are emitted:

Designation	Unit	Comment
Office BP		
Central systolic blood pressure (cSys)	mmHg	
Central systolic blood pressure (cDia)	mmHg	
Central pulse pressure (cPP)	mmHg	
Pulse pressure amplification		Not available in the USA
Hemodynamic		
Stroke volume (SV)	ml	Not available in the USA
Cardiac output (CO)	l/min	Not available in the USA
Total vascular resistance (TVR)	s·mmHg/ml or dyn·s/cm ⁵	Not available in the USA
Cardiac Index	l/min·1/m ²	Not available in the USA
Arterial Stiffness		
Augmentation pressure	mmHg	Not available in the USA
Augmentations Index@75 (Alx@75) [90 % CI*]	%	In the USA: In patients 40 years of age and older
Pulse wave velocity (PWV) [90 % CI*]	m/s	Not available in the USA

*Confidence Interval

Pulse wave analysis is performed in the office. Both individual and triple PWA measurements are available. Triple PWA measurement involves three consecutive PWA measurements with a brief pause between each of the measurements. This method is used to detect possible white coat hypertension in the doctor's office.

17.1 In-office pulse wave analysis

17.1.1 PWA license key for Mobil-O-Graph[®] / Tel-O-Graph[®] BT

A license key is required to perform pulse wave analysis using the Mobil-O-Graph[®] or Tel-O-Graph[®] BT. A license key is available from the manufacturer.

 **Note**

PWA measurement is only possible with the Bluetooth[®] interface.

Installation:**Pre-conditions:**

- The blood pressure monitor/measuring device and computer are switched on
- The Bluetooth[®] interface for the blood pressure monitor/measuring device is known by HMS CS (paired)
- The blood pressure monitor is connected with HMS CS via Bluetooth[®]

1. Once the Bluetooth[®] connection between HMS CS and the blood pressure monitor/measuring device has been established, click on the **Prepare device** button.

The **Prepare device** window opens.

2. Click on the **Activation** tab.

3. Enter the license key and click on **Send**.

A window stating that the PWA Flatrate license has been activated appears.

4. Click on **OK**, then on **Save**.

The PWA Flatrate license is now successfully activated.

 **Note**

The **Prepare device** window shows under the **Device test** tab whether PWA Flatrate is available for the device in question.

17.2 PWA license dongle for Mobil-O-Graph[®] (Not available in the USA)

The in-office PWA screening functionality is available in the Mobil-O-Graph[®] firmware from version 200007 onwards in conjunction with HMS CS version 2.0 or higher.

If you wish to use Mobil-O-Graph[®] to perform pulse wave analysis, you will need a USB license dongle. A license dongle is available from the manufacturer.

Installation:

1. Insert the PWA license dongle into a free USB socket on your computer.
2. Then start a PWA measurement.

 **Note**

24-hour PWA is not possible using the PWA license dongle.

17.3 In-office pulse wave analysis

17.3.1 Individual pulse wave analysis

 **Note**

Individual PWA measurements are possible using Mobil-O-Graph[®] (license key or dongle) / Tel-O-Graph[®] BT (license key).

Individual PWA measurement is only possible with the Bluetooth[®] interface.

Performing an individual pulse wave analysis:

Pre-conditions:

- The blood pressure monitor/measuring device and computer are switched on
- The interface for the blood pressure monitor/measuring device is known by HMS CS (paired)

1. Position the cuff of the blood pressure monitor/measuring device on the patient and connect the cuff to the device.

**WARNING**

Please note the warnings in the instructions for use of the blood pressure monitor/measuring device.

2. Select the relevant patient or sign-in a new patient in HMS CS.

The PWA measurement is always assigned to the patient currently open.

Note

To perform a PWA, the **age**, **size** and **weight** of the patient must be entered in HMS CS.

3. Connect the blood pressure monitor to HMS CS using **Bluetooth®**.
The **Device action** window opens.
4. Click on **PWA Measurement**.
The **PWA Measurement** window opens.
5. Click on **OK** to confirm that you have applied the blood pressure cuff.
The PWA measurement will start.
6. Once the entire PWA measurement process has been successfully completed, click on **Save**.

17.3.2 Triple pulse wave analysis

Note

Triple PWA measurement is possible using Mobil-O-Graph® (license key or dongle) / Tel-O-Graph®BT (license key).

Triple PWA measurement is only possible with the Bluetooth® interface.

Performing a triple pulse wave analysis:

Pre-conditions:

- The blood pressure monitor/measuring device and computer are switched on
- The Bluetooth® interface for the blood pressure monitor/measuring device is known by HMS CS (paired)

1. Position the cuff of the blood pressure monitor/measuring device on the patient and connect the cuff to the device.

WARNING

Please note the warnings in the instructions for use of the blood pressure monitor/measuring device.

2. Select the relevant patient or sign-in a new patient in HMS CS.

Triple PWA measurements are always assigned to the patient currently open.

Note

To perform a PWA, the **age**, **size** and **weight** of the patient must be entered in HMS CS.

3. Connect the blood pressure monitor/blood pressure measuring device to HMS CS using **Bluetooth®**.
The **Device action** window opens.

4. Click on Triple PWA measurement.
- The Triple PWA measurement window appears.
5. Start triple PWA measurement by clicking on **OK**. (This is to confirm that you have applied the blood pressure cuff.)

**Note**

If the battery charge level of the measurement device is insufficient for a measurement, the following warning symbol will appear:

**Note**

If the quality of the analysis is poor, an additional measurement will be automatically performed.

6. Once a triple PWA measurement has been performed, click on **Save**.
- A **PWA statistics** window appears, listing the three PWA measurements for comparison. The individual measurements are displayed in different colors.
7. Click on **Print** to print the statistics.
8. Close the **Statistics** window by clicking on **OK**.

The window is then no longer retrievable. The individual PWA measurements then appear consecutively in the list of measurements already taken.

17.4 24-hour PWA

Note

24-hour PWA measurement is possible with Mobil-O-Graph® and a PWA license key.

24-hour PWA measurement is only possible with the Bluetooth® interface.

17.4.1 Performing 24-hour PWA using Mobil-O-Graph®

Pre-conditions:

- The blood pressure monitor and computer are switched on
- The Bluetooth® interface for the blood pressure monitor/measuring device is known by HMS CS (paired)

1. Connect the blood pressure monitor to HMS CS using **Bluetooth®**.

The **Device action** window opens.

2. Once the Bluetooth® connection between HMS CS and the blood pressure monitor has been established, click on the **Prepare device** button.

The **Prepare device** window opens.

3. In day periods it is possible to set the number of PWA measurements.

4. Set the number of PWA measurements under the day periods. Options here include 30,15,12, etc., although these depend on the number of blood pressure measurements.

5. If the buttons in the **Prepare device** window are shown in red, click on them.

The corresponding settings are adjusted.

6. Save the settings by clicking on **Save**.

7. Position the cuff on the patient and connect the hose of the cuff to the blood pressure monitor.

WARNING

Please note the warnings in the instructions for use of the blood pressure monitor.

8. Ensure that the blood pressure monitor is working as desired, then press the START button

to start a manual measurement.

Note

Once the blood pressure measurement is complete, a PWA measurement will be performed with the device pumping up to diastolic pressure again and recording the pulse at that point.

17.4.2 Uploading and analysis of the 24-hour PWA

The data is uploaded and analyzed in the same way as conventional 24-hour ABPMs.

For more information, refer to the chapter "Transferring and storing long-term measurement results from the blood pressure monitor".

17.5 Displaying the pulse wave analysis

Once an in-office PWA has been successfully completed, the following analysis will be displayed:


Fig. 2: Representation of a pulse wave analysis in the HMS CS

Based on 10 pulse wave measurements, a filtered and average pulse wave is determined, which in turn is used to calculate the central aortic pulse wave.

The augmentation index (Alx) is repeatedly cited in medical literature as being dependent on gender, age and heart rate. This is why the use of a standardized presentation based on these parameters is preferred. A reference augmentation index is initially determined at a heart rate of 75 bpm using empirical regression⁶. This parameter is then known as Alx@75. If one were to examine a representative sample of the population as described in⁷, for example, the result would be an age-

⁶ Wilkinson I.B. et al. Heart Rate Dependency of Pulse Pressure Amplification and Arterial Stiffness. American Journal of Hypertension 2002; 15:24-30.

⁷ Fantin F. et al. Is augmentation index a good measure of vascular stiffness in the elderly? Age and Ageing 2007; 36: 43-48.

dependent estimate for Alx@75, plus an associated confidence interval. These pertinent studies have also revealed a significant difference in the average Alx@75 between men and women.

The average values displayed below with 90 % confidence intervals were determined on the basis of several internal studies⁸ with an internally collected representative sample of about 2000 people from the population. Similar to the previously mentioned studies, the measurements taken internally revealed an increase in the Alx up to the age of 55. Both genders then reached a plateau. The gender difference in Alx is between 8 and 10 %. Where measurement values exceed the gender-specific and age-specific interval, further studies in accordance with the European Treatment Guidelines for Hypertension⁹ are recommended in order to determine the reasons of the disorder.


Fig. 3: Mean value and 90% confidence interval for Alx@75

⁸ Nunan D. et al. Assessment of central haemodynamics from a brachial cuff in a community setting. BioMed Central 2012; 1471-2261

⁹ The Task Force for the Management of Arterial Hypertension of the European Society of Hypertension (ESH) and the European Society of Cardiology (ESC). 2007 Guidelines for the management of arterial hypertension. European Heart Journal 2007; 28: 1462-1536.

17.5.1 Alx and PWV history

To open the Alx and PWV history, click on the **Arterial stiffness** tab 


Fig. 4: The Alx history shows the Alx@75 depending on age.

18 Error messages

General

Error	Possible cause	Solution
Measurement data can no longer be opened/displayed in HMS CS.	An error occurred during patient data storage.	Delete the patient in question (toolbar) and recreate it.
No patient number	The device has not been initialized (i.e. the patient number was not transferred when preparing a measurement).	The patient number can be transmitted after the measurement. This will not adversely affect the measurement data.
There was a problem with the connection between the device and PC.	The wrong COM interface is set.	Set the correct interface.
Implausible data	This error message appears if the office software is attempting to open a patient and uses a wrong GDT ID in the process.	In this case, contact the customer support for the office software.
Wrong URL format	It is very likely that there is a typing error in the data source.	Check whether you have used \ instead of /.

Mobil-O-Graph®

Error	Possible cause	Solution
The time and date of the blood pressure monitor have not been updated and are different to that set in HMS CS.	<p>The blood pressure monitor has been stored for a long time without battery power.</p> <p>The internal backup battery of the blood pressure monitor is depleted. The metrological check of the blood pressure monitor that is required every two years has not been performed. (The backup battery is replaced during the metrological check.)</p>	<p>The date and time can be reset using HMS CS whenever the batteries are changed.</p> <p>Send the device for the metrological check to your specialist supplier or directly to IEM GmbH.</p>
There was a problem with the connection between the blood pressure monitor and PC.	<p>The wrong COM interface is set.</p> <p>Cable plug or socket is defective.</p>	<p>Set the correct interface.</p> <p>Inspect the plug and the socket on the blood pressure monitor. Ensure that the pins are not bent; bent pins may prevent contact.</p>
	The blood pressure monitor is not in transmission mode (the display shows the time).	Switch the blood pressure monitor off and then on again without removing the connection cable.
No measurements were conducted during the nocturnal phase.	<p>The batteries were prematurely depleted.</p> <p>The patient has switched off the blood pressure monitor.</p>	<p>The batteries may be defective (please contact your distributor).</p> <p>Draw the patient's attention to the need for a complete 24-hour measurement.</p>
No automatic measurements will be performed.	No manual measurement performed when fitted.	A valid measurement must always be performed manually when the device is fitted.

Tel-O-Graph[®]

Error	Possible cause	Solution
The time and date of the blood pressure monitor have not been updated and are different to the time and date set in HMS CS.	<p>The blood pressure monitor has been stored for a long time without battery power.</p> <p>The internal backup battery of the blood pressure monitor is depleted.</p> <p>The metrological check of the blood pressure monitor that is required every two years has not been performed. (The backup battery is replaced during the metrological check.)</p>	<p>The date and time can be reset using HMS CS whenever the batteries are changed.</p> <p>Send the device for the metrological check to your specialist supplier or directly to IEM GmbH.</p>

IEM[®]

HMS CS

Manual de instrucciones

Hypertension Management Software Client Server (HMS CS)

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³ Kazuaki Negishi et al. Importance of Calibration Method in Central Blood Pressure for Cardiac Structural Abnormalities. American Journal of Hypertension 2016.

⁴ Athanase D. Protopsaltis et al. Left-ventricular hypertrophy is associated better with 24-h aortic pressure than 24-h brachial pressure in hypertensive patients: the SAFAR study. American Journal of Hypertension 2014.

⁵ Thomas Weber et al. Validation of a Brachial Cuff-Based Method for Estimating Central Systolic Blood Pressure. Journal of the American Heart Associations 2011.

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