

## Declaration of Conformity

Wir erklären in eigener Verantwortung,  
dass das Produkt:

Mobil-O-Graph NG®  
Norav Branding: NBP-24 NG

UMDNS: 16-173  
Elektronischer, nicht-invasiver Blutdruckmonitor  
die Anforderungen der Richtlinien erfüllt:

*Richtlinie 93/42/EWG (MDD)*  
*Richtlinie 2014/53/EU (RED)*  
*Richtlinie 2011/65/EU (RoHS)*

### Konformitätsbewertungsverfahren

nach Anhang II.3, Richtlinie 93/42/EWG  
incl. Änderung durch die Richtlinie 2007/47/EG

### Klasse IIa

nach Anhang IX, Regel 10, Richtlinie 93/42/EWG

### Benannte Stelle (nur RL 93/42/EWG)

TÜV NORD CERT GmbH  
Langemarckstr. 20  
45141 Essen, Deutschland  
Kennnummer: CE 0044

### Referenznormen

Alle anwendbaren Prüfungen nach  
IEC 60601-1:2005+A1:2012, IEC 60601-1-2:2014,  
IEC 60601-1-11:2015, IEC 80601-2-30:2018  
sind erfolgreich durchgeführt worden.  
Die komplette Liste der anwendbaren Normen ist Bestandteil der  
Technischen Dokumentationen und kann  
auf Wunsch eingesehen werden.


### Harmonisierte Normen RED

EN 300 328 v2.2.2

### Gültig bis: 2024-05-24

basierend auf dem EG-Zertifikat  
04 232 041526 ed. 5

Place and date: Aachen, 2021-04-01



Thomas Korth (CEO)

We hereby declare in our own responsibility  
that the device:

Mobil-O-Graph NG®  
Norav Branding: NBP-24 NG

UMDNS: 16-173  
electronic, non-invasive blood pressure monitor  
is in compliance with the following directives:

*Directive 93/42/EEC (MDD)*  
*Directive 2014/53/EU (RED)*  
*Directive 2011/65/EU (RoHS)*

### Conformity Assessment Procedure

according to *annex II.3* of Council Directive 93/42/EEC  
as amended through Directive 2007/47/EC

### Class IIa

according to Annex IX, Rule 10, Directive 93/42/EEC

### Notified Body (93/42/EEC only)

TÜV NORD CERT GmbH  
Langemarckstr. 20  
45141 Essen, Germany  
Identification number: CE 0044

### Reference Standards

All applicable tests according to  
IEC 60601-1:2005+A1:2012, IEC 60601-1-2:2014,  
IEC 60601-1-11:2015, IEC 80601-2-30:2018  
have successfully been performed.

The full list of applicable standards is part of the  
Technical Documentation and can be provided  
upon request.

### Harmonized Standards RED

EN 300 328 v2.2.2

### Valid until: 2024-05-24

based on the EC Certificate:  
04 232 041526 ed. 5

